
# usmile

## User Manual

## usmile Sonic Electric

## Toothbrush Y10

Get more when you  
order usmile!

Follow us @usmile

TikTok

@usmile.official

@usmile.global

# Contents

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# Environmental Protection

All products bearing this symbol are waste electrical and electronic equipment (WEEE as in directive 2012/19/EU) which should not be mixed with unsorted household waste.

Instead, you should protect human health and the environment by handing over your waste equipment to a designated collection point for the recycling of waste electrical and electronic equipment, appointed by the government or local authorities. Correct disposal and recycling will help prevent potential negative consequences to the environment and human health.

Please contact the installer or local authorities for more information about the location as well as terms and conditions of such collection points.

Warning: Only remove the rechargeable battery when you discard the appliance. Make sure the battery is completely empty when you remove it.

This product is charged by Type-C charging cable, adapter output voltage/current is 5Vdc/1A. If you use a third-party charger, the recommended output voltage/current of the adaptor is 5Vdc/1A, and the adapter shall be CE approval type. Using unauthorized charger may cause danger.

Hereby, our company of Guangzhou Stars Pulse Co., Ltd., declares that this product is in compliance with the applicable Directives and European Norms and amendments.

This product is charged by Type-C charging cable, adapter output voltage/current is 5Vdc/1A. If you use a third-party charger, the recommended output voltage/current of the adaptor is 5Vdc/1A, and the adapter shall be UKCA approval type. Using unauthorized charger may cause danger.

Hereby, our company of Guangzhou Stars Pulse Co., Ltd., declares that this product is in compliance with the applicable Directives and United Kingdom Norms and amendments.

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Výrobce: Guangzhou Stars Pulse Co., Ltd.